

# **NUVERIA LABS GELLING FIBER DRESSINGS**

Gentle care for a wide range of wounds.





#### **DESIGNED FOR**

# **Wound Management:**

- Wide range of partial to full-thickness wounds such as pressure injuries, foot and leg ulcers, lacerations and abrasions
- Wounds with moderate to heavy drainage
- Reinforced ribbon dressing can be used in deep and tunneling wounds

### **BENEFITS**

This dressing is comprised of carboxymethylcellulose (CMC) fibers, which transform into a gel when they come in contact with fluid or wound drainage.



#### **ABSORBENT**

Highly absorbent\* and able to keep drainage away from the wound.



#### **STRONG**

Structurally solid while being highly conformable.



#### **GENTLE**

Promotes the removal of nonviable tissue from the wound. while maintaining healthy tissue. Allows for gentle application and atraumatic removal.



## REMEDIAL

Creates and maintains an optimal moist wound healing environment.







# Gelling fiber dressings SKU details:





Size	Qty per box	Qty per case
2" x 2" (5x5cm)	10	100
4" x 5" (10x12.7cm)	10	100
6" x 6" (15x15cm)	5	50
1" x 18" (2.5x45.7cm)	5	50
	2" x 2" (5x5cm)  4" x 5" (10x12.7cm)  6" x 6" (15x15cm)	2" x 2" (5x5cm) 10  4" x 5" (10x12.7cm) 10  6" x 6" (15x15cm) 5

For detailed product information, including indications for use, precautions, and warnings, please consult the product's IFU.



# CLICK OR SCAN TO ORDER TODAY

Fast, free delivery with Amazon Prime.

No minimum orders. Business pricing available based on monthly usage.